

## Advisories

- Use a dry cloth or alcohol free wet wipes to clean. Do not rinse under water as the mask is not waterproof. DO NOT use any other cleansing products other than those stated. Use of other cleansing products may damage the +maskology **PHOTON LED MASK**.
- Store the device in a cool, dry space out of direct sunlight or high temperatures.
- Keep out of reach of children.
- Never attempt to open or repair the device yourself.
- The straps contain a mixture of polyester, nylon and lycra. If you are allergic to any of these materials, you may not be able to use the straps.
- The mask cannot be used while it is in the charging mode.
- Store the +maskology **PHOTON LED MASK** flat and at room temperature in dry conditions, away from direct sunlight.

## Advice/Warnings

- DO NOT use it without first removing creams, serums, lotions, and make-up BEFORE treatment. Some ingredients in these products may cause light sensitisation or may reduce the effect of the treatment. If you use The +maskology **PHOTON LED MASK** without first removing these products, you may cause a skin reaction, damage your device, or you may reduce the effect of the LED treatment.
- DO NOT attempt to modify The +maskology **PHOTON LED MASK** in any way. Modifying The +maskology **PHOTON LED MASK** will void your warranty. Modifying The +maskology **PHOTON LED MASK** could cause an electric shock when the device is used again. The device contains no serviceable parts.
- TURN OFF and REMOVE The +maskology **PHOTON LED MASK** immediately if it feels hot or if your skin feels uncomfortably warm. The The +maskology **PHOTON LED MASK** does not produce any heat.
- DO NOT use the +maskology **PHOTON LED MASK**
  - If you are particularly sensitive to light.
  - If you are pregnant, planning on becoming pregnant or breastfeeding. The risk to the foetus, neonate or pregnant women is unknown
  - If you suffer from epilepsy
  - If you are currently using medication that can cause skin sensitivity to light
  - When sleeping, walking or driving.
  - If you suffer from light induced headaches.
  - To treat any other conditions apart from those listed in the indications for use.
  - If you suffer from lupus erythematosus, photosensitive eczema or albinism. If you use it to treat lupus erythematosus, photosensitive eczema or albinism you may cause a severe skin reaction.
  - If you suffer from any photosensitive disorder (sensitisation to light). If you use the mask and you suffer from a photosensitive disorder you may cause a severe skin reaction.
  - If you are taking any medication that can cause photosensitivity. If you use the mask and you are taking any medication that can cause photosensitivity you may cause a severe skin reaction.
  - If you suffer from any genetic conditions of the eye. Photosensitivity is a common side effect of various medications. These can include certain antibiotics, chemotherapy drugs, and diuretics. If you are unsure about any medication you may be taking, consult your healthcare provider.
  - In any way that is not exactly as specified in these instructions for use, other than exactly as specified in these instructions for use has not been tested and the risk is unknown.
  - If you think it may have come into contact with water. Keep away from HOT surfaces.
  - Near water or anywhere water can be spilt onto the equipment.
  - If the device, battery, or charging plug are damaged in any way.